

VITEK 2™ – compact 15



Automated Solutions

With a 15-card capacity, the VITEK® 2 Compact 15 offers the same design as the VITEK 2 Compact 30 and the VITEK 2 Compact 60. VITEK brings proven automation and recognized reliability to every laboratory.



BIOMÉRIEUX

Each of the VITEK® 2 family of instruments is based on the same innovative VITEK 2 technology that includes: an expanded identification database, the most automated platform available, rapid results, improved confidence, and limited training time.



Bringing VITEK 2 productivity and automation to every laboratory

- **Proven VITEK 2 Technology** available for every size laboratory
- **Proven identifications** with advanced colormetric optics that allow for an extended database and increased online identification
- **Easy training and implementation** reduce training time with icon-driven, user-friendly Windows-based PC software
- **Time, labor, and money savings** improve workflow through automation and barcodes
- **Peace-of-mind** and increased confidence with AES, the only second-generation expert system available for antibiotic susceptibility testing. AES provides the added value of a red, yellow or green light and improves time management for critical results
- **Same capabilities in all labs** expand the possibility for in-house patient care

Specifications

Weight:

VITEK 2 Compact 15/30/60
(15/30/60 card capacity) 75 kg (165 lb)

Electrical Power Requirements:

- 100/120 VAC (50-60 HZ)
- 220/240 VAC (50-60 HZ)

Heat Dissipated:

VITEK 2 Compact: 1025 BTU/Hr. (nominal)

Environmental Requirements:

- Operating Ambient Temperature Range: 59°F-86°F (15°C to 30°C)
- Operating Humidity Range: 20% to 80% relative humidity, non-condensing

Altitude:

Up to 2000 m.

Capacity:

15, 30 or 60 cards (2 modules can be connected to the same computer).

Instrument Dimensions:

H 23.6"
W 28.6"
D 26.8"

Contact your local bioMérieux representative to find out more about **VITEK 2 Compact 15**.

bioMérieux, Inc.
100 Rodolphe Street
Durham, NC 27712
U.S.A.
Tel: (800) 682 2666
Fax: (800) 968 9494

bioMérieux Canada, Inc.
7815 Henri-Bourassa Boulevard West
St. Laurent (Québec)
Canada H4S 1P7
Tel: (1) 514 336 7321
Fax: (1) 514 807 0015

www.biomerieux-usa.com



B I O M É R I E U X