

Complete solution for COVID-19 response


New standard for COVID-19 test to overcome the evolving SARS-CoV-2

[Learn More](#)

[▶ Video](#)

Primary screening

for SARS-CoV-2 & Flu and RSV in one tube




Allplex™
SARS-CoV-2/FluA/FluB/RSV Assay

[PDF](#)

[▶ Video](#)

Selective combination test

for SARS-CoV-2 & common respiratory viruses



Allplex™
SARS-CoV-2 Assay

Allplex™
RV Essential Assay

[PDF](#)

[▶ Video](#)

Features

	Allplex™ SARS-CoV-2/FluA/FluB/RSV Assay	Allplex™ SARS-CoV-2 Assay
No. of target genes for SARS-CoV-2	3	4
Syndromic test	One tube (Flu A, Flu B, RSV A/B)	Selective combination test
Automated platform	✓	✓
Extraction-free PCR	-	✓
Self-collected sample	✓	-
Test process verification	Endogenous & Exogenous dual IC	Exogenous IC